

# Hybridizing sample for 60k

## Prepare the 10× Blocking Agent:

1. Add 1,350  $\mu\text{L}$  of DNase/RNase-free distilled water to the vial containing lyophilized 10× aCGH Blocking Agent (included in the Oligo aCGH/ChIP-on-chip Hybridization Kit).
2. Leave at room temperature for more than 6 hours and mix on a vortex mixer to reconstitute sample before use or storage.

## Hybridizing sample

1. Check if the incubator is at 95°C
2. Re-hydrate pellets with 11.9 $\mu\text{L}$  of ddH<sub>2</sub>O (put on the wall of the tube).
3. Spin down, keep in dark.
4. Prepare Hybridization Solution Master Mix

60k	Each	1 slide	2 slides	3 slides	4 slides
2*HI-RPM Hybridization Buffer	27.5	247.5	495	742.5	962.5
10*Acgh Blocking Agent	5.5	49.5	99	148.5	192.5
Formamide	5.5	49.5	99	148.5	192.5
Cot-1 DNA	2.4	21.6	43.2	64.8	84
Universal standard	2.2	19.8	39.6	59.4	77
total buffer volume	43.1	387.9	775.8	1163.7	1508.5

\* Add US in dark room

5. Add 43.1  $\mu\text{L}$  of Hybridization Solution Master Mix into the tube
  6. Mix well (~15 seconds) and spin down
  7. Denature labeled DNA at 95°C for 3 minutes
  8. Immediately transfer tubes to a 37°C incubator. Incubate at 37°C for 30 minutes
  9. Spin 1 minute at 6000  $\times$  g, keep in 37°C incubator
  10. Hybridization
- Check if hybridization oven is setting at 67°C and 20 rpm.

- Load a gasket slide into the Agilent SureHyb chamber base with label “Agilent” facing up

- Take 50  $\mu$ L of the mixture and load about 48  $\mu$ L onto a gasket well

**Caution: Do not touch the gasket slides when loading the hybridization sample mixture.**

- Put a microarray slide with label “Agilent” down onto the gasket slide, the numeric barcode side is facing up.

- Put the SureHyb chamber cover onto the sandwiched slides.

- Slide the clamp assembly onto both pieces.

- Hand-tighten the clamp firmly onto the chamber.

- Vertically rotate the assembled chamber, make sure that bubbles can smoothly move.

- Tap the assembly on a hard surface if necessary to move stationary bubbles.

- Place assembled slide chamber in the hybridization oven.

**Caution: Be sure to balance the loaded hybridization chambers on the rack.**

11. Hybridize the sample for 22-24 hours at 67°C.

12. Pre-warm 200mL **Wash Buffer 2** at 43°C overnight.

### Washing Arrays

	Dish	Wash buffer	Temperature	Stirrer	Time
Disassembly	#1	Wash Buffer 1	RT		
1st wash	#2	Wash Buffer 1	RT	200 rmp	5 min
2nd wash	#3	Wash Buffer 2	37 °C	140 rmp	1 min

1. Minimize exposure of the slide to air.

2. Touch only the barcode portion of the microarray slide or its edges.

3. Spin dry the slide.

### Scanning

1. Insert each slide with label “Agilent” side up and barcode end first into the slide holder.

2. Scan using the Multi-TIFF settings.