

Labeling of microbial community DNA for microarray hybridization

1. Prepare dNTP mix: (using 100 mM dNTPs)

5 μ l each of dA/G/CTP

2.5 μ l of dTTP

82.5 μ l DEPC Water

2. Prepare the DNA/ Random primer mixture:

Random primer (Life Technologies, random hexamers, 3 μ g/ μ L) 5.5 μ L

Target or control DNA 29.5 μ L

Total volume 35 μ L

(1) Transfer 5.5 μ L of Random primer to a PCR tube.

(2) Transfer 29.5 μ L DNA to the Random primer. Usually, 500 ng to 1000 ng of DNA should be used. Adjust volume to 29.5 μ L with water; or if the original concentration of a sample is very low, the DNA can be concentrated.

(3) Mix DNA and Random primer thoroughly and incubate at 99.9 $^{\circ}$ C for 5 min (using thermocycler)

(4) Immediately chill tubes on ice.

3. Prepare labeling premix

In a separate tube, for each reaction combine:

10X buffer (included with the klenow)	5 μ L
dNTP mix	2.5 μ L
Klenow (imer)	1 μ L
CyDye* (25 nM)	0.5 μ L
H ₂ O	6 μ L
<hr/>	
Total	15 μ L

* CyDyes are light sensitive – dispense in the dark room. This can be added last

(1) Transfer 15 μ L of the labeling premix to the DNA/ Random primer mixture and mix well.

(2) Incubate the reaction at 37 $^{\circ}$ C for 6 hrs.

(3) Heat-inactivate the enzyme by incubating the reaction at 95 $^{\circ}$ C for 3 min and then cool at 4 $^{\circ}$ C

*Labeled DNA can be kept at 4 $^{\circ}$ C in the thermocycler if labeling is done overnight.

4. Purification

- (1) Purify labeled DNA with Qiagen QIAquick Kit as specified by the manufacturer;
- (2) Elute the DNA using 100 μL H_2O or EB buffer
- (3) Check CyDye incorporation using Nanodrop

Minimum dye incorporation: $\text{pmol} > 50$ ($\text{pmol}/\mu\text{L} * \text{total } \mu\text{L}$)

[for example: Labeled DNA eluted with 100 μL of EB buffer and the $\text{pmol}/\mu\text{L}$ is 0.8; $100*0.8=80$ – labeling is fine, proceed with hybridization]

5. Dry using a Speed Vac

Program: Temp. 45 $^{\circ}\text{C}$, Run 2 hours, Vacuum (Level) 5.1 (preset)

Press “Auto Run”

Pellet should be dry at the end of the run.